

UPGRADE YOUR MEAT



Enhanced tenderness

for a given cooking regime and original quality of meat

Add 0.4% - 1.4% baking soda by weight of meat
to a 30 minute marinade prior to cooking

Bicarbonate soaking results in higher pH and water
binding rendering the connective tissue more susceptible
to heat denaturation resulting in more tender beef



**Arm & Hammer® Baking Soda
1 lbs. Box**

Item #33200 - 84104
Pack: 24/16 oz
UPS Shippable



**Arm & Hammer® Baking Soda
2 lbs. Box**

Item #33200 - 01140
Pack: 12/32 oz
UPS Shippable



**Arm & Hammer® Baking Soda
4 lbs. Box**

Item #33200 - 01170
Pack: 6/4 lbs
UPS Shippable



**Arm & Hammer® Baking Soda
13.5 lbs. Resealable Pouch**

Item #33200 - 01961
Pack: 80/13.5 lbs